

Monopolar Radiofrequency

Fight against Tissue Laxity



Improves your face and body's skin

The monopolar radio frequency is a device that guarantee the skin rejuvenation and fight against the skin laxity of face and décolletè skin, and many others body part like abdomen, hips, gluteus and triceps.



DATA SHEET

Descrizione	Caratteristiche
Alimentazione	220 - 240Vac 50/60Hz
Assorbimento	50 VA
Fusibili	2 x T 1 A
Potenza max RF	50 W
Frequenza Lavoro	433 KHz
Visualizzazione parametri	Tramite display grafico
Impostazione dei parametri	Guidata da menu ad icone e grafica
Comandi	5 Pulsanti : up – down - dx - sx - conferma
Allarmi ed indicazione	Temperatura dissipatori elevata I
Sicurezze	Circuiti di retroazione Alimentatore protetto contro i cc
Controlli	Controllo a microprocessore di tutti i parametri e del funzionamento dell'alimentatore
Temperatura di funzionamento	Da 10 a 40°C
Umidità rel. Di funzionamento	Da 30 a 75%
Tipo Apparecchio	Tipo BF
Classe Apparecchio	Classe I
Grado di protezione	IPX0


FACE'S TREATMENT

The treatment is recommended also for the face's skin, the **results** are evident from the **very first treatments**.



BODY'S TREATMENT

Working on the body, **the results of the radio frequency are evident and stable over time**. The device improves tissue drainage and increase the microcirculation of the tissue helping with the elimination of the first stages of cellulitis.





RE-CO s.r.l.

ITALY : Via Trento, 16 Canonica d'Adda (BG) 24040

Tel. [02/90989319](tel:0290989319) Fax. 02/90988626 Mail info@recoweb.it