

# Air Cooling System

**New revolution!!!** The most powerful Air Cooling System



Model : Air Cooling System

- **Convenient and economical for users**
  - It is very economical because treatment times can be considerably reduced.
  - It is very convenient because using is very easy
  - It is very useful because consumable parts are not needed.
- **Very effective and safe for patients**
  - It can reduce skin damage because it can be used with lasers.
  - It can reduce pain and recovery time.
  - It can protect skin because cooling can be made evenly.
- **The suitable lasers**
  - Fractional CO<sub>2</sub> Laser, IPL, RF Microneedle type, Er yag, Diode laser, Alexandrite, Nd:YAG(Q-Switched, long pulsed yag)
- **Special Features**
  - Air flow controlling system
  - Cooling air temperature (Max. - 35 degree)
  - Quick cooling with high quality compressor
  - Auto adjustable controlling system
  - Auto detecting system to keep best condition



## TECHNICAL SPECIFICATIONS

A.C Input	220~240V, 50/60Hz
Air Level	10 Step
Size	400(W)×610(L)×963(H)
Air Flow Output	1000 L/Min
Temperature	0~-35°C
Weight	50Kg

