


Dermoabrasion

Fights the imperfections of age



The dermabrasion helps in “smooth” the edges of the scars, independently from the origin (post-acne, post-trauma, post-surgery)

The treatment with the dermabrasion is suggested in various circumstances. It can be applied both face and body.



ELIMINATES

↳ Acne and acne scars.

↳ Wrinkles.

↳ Expression lines.

↳ Dark circles and bags

↳ Dull Skin, thick, with unhealthy complexion.

↳ Sunspots.

↳ Various discoloration, hormonal, post-partum, solar, age.

↳ Stretch Marks.

↳ Scars.

↳ Tatto to eliminate



This treatment is carried out with special machines equipped with a particular handpiece that is placed on the part of the face or body to be treated.

Solid particles come out of this handpiece at high speed and have the task of removing the superficial layers of the skin, exfoliating it.

The blemishes that are best suited to be corrected with the dermabrasion technique are:

- Depressed post-acne scars of the face
- Post-surgical scars of the face
- Post-traumatic scars of the face
- Static and expression wrinkles of the upper lip
- Seborrheic and solar keratoses



DATA SHEET



Descrizione	Caratteristiche
Alimentazione	230Vac 50/60Hz
Assorbimento	70 VA
Fusibili	2 x T 1,5 A
Potenza max Vacuum	Regolabile da 0 a -800 mBar
Visualizzazione parametri	Tramite Vacuometro meccanico
Impostazione dei parametri	Con manopola
Comandi	manuale
Sicurezze	Fusibili
Temperatura di funzionamento	Da 10 a 40°C
Umidità rel. Di funzionamento	Da 30 a 75%



RE-CO s.r.l.

ITALY : Via Trento, 16 Canonica d'Adda (BG) 24040

Tel. [02/90989319](tel:0290989319) Fax. 02/90988626 Mail info@recoweb.it