

DECOBLEND

Professional Blend Electrolysis technology

Made in Italy technology for hair removal



Permanent Hair Removal Beyond the Limits of Laser

DECOBLEND is Decomedical's professional Blend electrolysis equipment, developed to provide reliable, progressive and permanent hair removal through the automatic combination of radiofrequency thermolysis and low-voltage direct current electrolysis.

Unlike light-based technologies, DECOBLEND effectively treats all hair colours, including white, blonde, grey and red hair, making it the ideal solution for completing every professional hair removal service.

Designed and manufactured in Italy, it combines proven Blend technology, treatment precision and long-term reliability for beauty professionals.

Why choose Decoblend

- Permanent hair removal
- Suitable for every hair colour
- Suitable for all skin phototypes
- Facial and body treatments
- Ideal around tattooed areas
- Made in Italy



Accessories



Two Technologies. One Permanent Solution.

DECOBLEND combines two complementary treatment methods within the same treatment cycle.

Thermolysis

Radiofrequency energy generates controlled heat inside the follicle, weakening the hair structure through a precise thermal action.

Electrolysis

A low-voltage direct current creates a chemical reaction that progressively destroys the follicular structures responsible for hair growth.

Blend Technology

The simultaneous action of both energies combines the advantages of each method, delivering highly effective and long-lasting permanent hair removal with exceptional treatment precision.

Why Professionals Choose DECOBLEND

- Permanent hair removal regardless of hair colour
- Effective on white, blonde, grey and red hair
- Suitable for every skin phototype
- Safe treatment around tattooed areas
- High-precision individual hair treatment
- Progressive and long-lasting results
- Comfortable treatments with controlled energy delivery

TREATMENT APPLICATIONS

Designed for Precision

DECOBLEND is suitable for progressive permanent hair removal on every facial and body area.

FACE

- Upper Lip
- Chin
- Eyebrows
- Sideburns
- Neck

BODY

- Underarms
- Bikini Area
- Breast
- Abdomen
- Arms
- Legs
- Fingers
- Toes

SPECIAL APPLICATIONS

- White Hair
- Blonde Hair
- Grey Hair
- Red Hair
- Hair Around Tattoos



**BEFORE AND AFTER
EXAMPLE**

HOW THE TREATMENT WORKS


Precise. Controlled. Professional.

Using the ergonomic needle-holder handpiece, a sterile disposable needle is gently inserted into the natural hair follicle.

Energy delivery is activated through the foot pedal, allowing complete operator control throughout the treatment.

At the end of the programmed cycle, the treated hair is gently removed while preserving the surrounding tissue.

The combination of precision and controlled energy delivery ensures reliable, progressive results with maximum treatment accuracy



BUSINESS OPPORTUNITY

Complete Every Hair Removal Portfolio

Laser technology has transformed professional hair removal, but it cannot effectively treat hairs with little or no melanin.

DECOBLEND fills this gap by enabling professionals to permanently remove white, blonde, grey and red hair while safely treating delicate and tattooed areas.

A Valuable Opportunity

- Expand your treatment portfolio
- Complete laser treatment plans
- Attract new clients
- Increase repeat treatments
- Differentiate your clinic
- Become the White Hair Removal Specialist

